STORZ

KARL STORZ—ENDOSKOPE

Endo-Arthroflator VET 62432520

Endo-Arthroflator VET 62432520

Endo-Arthroflator VET 62432520

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Instruction manual

Manual de instrucciones

Endo-Arthroflator VET 62432520

Endo-Arthroflator VET 62432520

Endo-Arthroflator VET 62432520

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.


This instruction manual is intended as an aid to correct installation, connection and operation of the CO₂ Endo-Arthroflator VET. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment.

Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

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Controls, displays, connectors, and their uses

- ① Main power switch
- ② Book symbol. Explanation: Please read the information provided with the product (e.g. the instruction manual)
- ③ Pressure status
 - The symbol is illuminated red when the input pressure is too low (< 10 bar) (with hysteresis)
 - The symbol is illuminated green when the input pressure is within the permissible range (≥ 10 bar)
- ④ Digital indication of patient-line setpoint/current actual pressure (setpoints are indicated by a dot to the right of the numbers)
- ⑤ Digital indication of setpoint/current actual gas flow rate (setpoints are indicated by a dot to the right of the numbers)
- ⑥ Digital gas volume consumption indicator
- ⑦ Insufflation line fitting, line to patient
- ⑧ RESET button, resets gas volume consumption indicator
- ⑨ Insufflation START/STOP button
- ⑩ Button to enter the gas flow setpoint. The fixed values are gradually increased by pressing the button (up to a maximum value, the minimum value then reappears when pressed again)
- ⑪ \pm Buttons for entering the actual patient-line pressure
- ⑫ ARTHRO mode symbol. Is illuminated when the ARTHRO mode is activated
- ⑬ Connection for gas supply (CO_2 bottle)
- ⑭ Potential equalization connector
- ⑮ Line fuses: 2 x T1AL/250 V
- ⑯ Power cord socket

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 Read the instructions carefully before operating the equipment!

 Gas input pressure display

 Patient pressure

 Insufflation pressure

 Gas Flow

 CO₂ Volume

 Gas outlet to the patient

 Gas volume reset

 START/STOP

 Potential equalization connector

 Alternating current

 Electronic information product pollution control (China RoHS)

 **DANGER:** Risk of explosion if used in the presence of flammable anesthetics.

 **CAUTION:** To reduce the risk of electrical shock, do not remove cover. Refer servicing to qualified personnel.

Keep out of reach of patients. Do not store liquids on or above the unit.

Refer to instruction manual for disinfection/sterilization instructions.

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Photographs of the unit	5
Controls, displays, connectors, and their uses	6
Symbols employed	7

General information

Description of the unit	10
-------------------------------	----

Safety instructions

Warnings and cautions	11
Intended use	14
Indications for use/Contraindications for use	14
User qualification	16
Safety precautions at the site of installation	17
Safety precautions when operating the unit	17
Safety features, acoustic signals	18

Installation and operating instructions

Unpacking the equipment	21
Basic equipment	21
Installing and connecting up the unit	21
Connection to a CO ₂ gas bottle.	23
Operating	24
Switching on LAP/ARTHRO mode	24
Preparing for CO ₂ insufflation	26
Initiating CO ₂ insufflation	28
Service functions.	29

Maintenance

Fuse replacement	30
Cleaning, disinfection and sterilization	31
Maintenance	36
Servicing and repair	36
Repair program	37
Limitation of liability	38
Manufacturer's warranty	38

Technical description

Troubleshooting	39
Technical data	40
Technical documentation	43
Plan of electronics	44
Plan of pneumatics	45

Spare parts, recommended accessories

List of spare parts	46
Accessories	47

Appendix

Cleaning agents and disinfectants	48
Electromagnetic Compatibility (EMC) Information	50

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Description of the unit

The Endo-Arthroflator VET 62432520 is an insufflator which is intended exclusively for use in veterinary medicine.

The Endo-Arthroflator VET has two operating modes: Both laparoscopic and arthroscopic operations can be performed with this device.

Switching on LAP/ARTHRO mode: When switched on, the device works exclusively in LAP mode. To access the ARTHRO mode, the user must deliberately switch over the device (actuate the power switch ① and, at the same time, press and hold down the + and - pressure button ⑪ for at least 3 seconds - the ARTHRO display appears). Confirm the flashing setpoint display by pressing one of the buttons ⑩ or ⑨.

The gas is supplied by a CO₂ gas bottle of medical quality.

Warning: Gas bottles may only be operated when upright. The connection must be at the top to supply the gas.

Warning: Consumption of large quantities of gas can cause the pressure reducer (in the device) to become iced which, in turn, may result in the device malfunctioning.

The patient pressure in LAP mode can be set in increments of 1 mmHg between 0 and 15 mmHg and in the ARTHRO mode in increments of 5 mmHg between 50 and 100 mmHg.

In the LAP mode, the flow can be set to 0 / 1.5 / 5 / 10 and 20 l/min. In the ARTHRO mode, the following values can be set 0 / 1.5 / 3 / 5 and 10 l/min.

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Warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words, **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

Warning: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

Caution: A Caution indicates that particular procedures or precautions must be followed to avoid possible damage to the product.

Note: A Note indicates special information about operating the product, or clarifies important information.

Warning: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

Note: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

Warning: The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.

Warning: Keep out of reach of patients.

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Warning: The Endo-Arthroflator VET is only suitable for use in veterinary medicine. It may not be used in human medicine.

Warning: The instructions and interface specifications for medical devices (e.g. CO₂ gas filter) and/or system components used in combination must be observed precisely.

Warning: Combinations of medical devices are only assured to be safe if they are identified as such in the respective instruction manuals or the intended purpose and interface specifications of the devices used in combination permit this (cp. IEC 60601-1-1).

Warning: Grounding reliability can only be achieved when the equipment is connected to "Hospital Only" or "Hospital Grade" outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

Warning: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

Warning: DO NOT use the unit in the presence of flammable anesthetics.

Warning: Always unplug the unit before all maintenance work (e. g. cleaning).

Warning: To avoid the risk of electrical shock, do not open the unit. Refer servicing to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

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Caution: Before connecting the unit to the electrical supply lines, verify that the line voltage stated on the unit's identification plate agrees with that of the electrical supply lines to be used.

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Caution: Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

Caution: To avoid contamination of the unit as a result of the reverse flow of gas or body fluid, a sterile filter (031122-25) must be inserted between the insufflation line and unit connection.

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Caution: If, during an arthroscopic intervention, an irrigation tube is also connected to the arthroscopy sheath, it must be ensured that the stopcock for the CO₂ supply on the arthroscopy sheath is closed before opening the stopcock for the irrigation fluid, in order to prevent a reverse flow of irrigation fluid into the unit which would damage it.

Warning: Electric installation at the installation site must conform with applicable national regulations. A connector with protective earth is essential in order to use the Endo-Arthroflator VET.

Warning: The unit may only be used indoors in dry areas. Do not use outdoors, in rooms with openings to the outside that cannot be closed off or in wet areas (bathrooms, etc).

Warning: Do not operate the unit in the presence of flammable gases or dust. Risk of explosion!

Warning: Always disconnect the unit prior to any maintenance work.

Warning: Danger of electric shock! Do not open. Refer repairs to authorized personnel only.

Caution: The Endo-Arthroflator VET may only be operated with tubing sets and accessories which have been designated as suitable for the device by KARL STORZ.

Caution: Care must be taken that liquids do not penetrate the housing. Do not store liquids on or above the unit.

Caution: Operate the unit with the voltage indicated on the provided label.

Caution: Use only fuses of the correct rating.

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Note: Make certain to observe national regulations/laws regarding disposal.

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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Intended use

The Endo-Arthroflator VET is a CO₂ insufflation device which has been exclusively approved for use in veterinary medicine. It is used for laparoscopic and arthroscopic examinations and operations. Only CO₂ gas which the manufacturer has designated as suitable for medical purposes may be used.

Use of the unit in fields other than those indicated above is not allowed for safety reasons.

The Endo-Arthroflator VET may only be operated with tubing sets and accessories which have been designated as suitable for the device by KARL STORZ.

Caution: If tube systems other than those prescribed are used, no guarantee can be given for safe operation of the instrument.

Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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Indications for Use/ Contraindications for Use Indications for Use

1. The Endo-Arthroflator VET provides CO₂ gas distension of the abdomen for diagnostic and/or operative laparoscopy.
2. The Endo-Arthroflator VET provides CO₂ gas distension of the abdomen during diagnostic and/or operative arthroscopy.

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Contraindications for Use

Use of the Endo-Arthroflator VET is contraindicated whenever laparoscopy or arthroscopy is contraindicated.

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Metabolic Acidosis and Resultant Cardiac Irregularity.

Warning: Prolonged intraabdominal pressures greater than 20 mmHg should be avoided. This can cause any of the following:

- Decreased respiration with compromised diaphragmatic excursion
- Decreased venous return
- Decreased cardiac output
- Acidosis

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High intra-abdominal pressure

Excessive absorption of CO₂ results from either excessive gas flow and/or excessive gas pressure. The abdomen can be adequately distended by the maximum pressure of 15 mmHg, which is valid in LAP mode. Little intravasation should occur at these levels.

Adequate respirations help avoid problems related to CO₂.

Note: At a pressure of 55-100 mmHg (LAP mode) or 120-180 mmHg (ARTHRO mode), the mechanical safety valves provide for pressure relief, thereby eliminating the risk of endangering the health of the patient.

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- **Hypothermia.** High flow rate insufflators without a gas flow heater and humidifier may present a potential risk for hypothermia.

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Precautions

To prevent patient cross contamination, a bacterial hydrophobic filter must be used.

Note: This filter (Art. No. 031122-25, mtp, for single use) also serves to prevent contamination and stops the device from being filled with irrigation liquid!

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
User qualification

The Endo-Arthroflator VET may only be used by vets and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the device.

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Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones.

This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using it (see the section "Testing the unit for proper operation"). During treatment using the Endo-Arthroflator VET the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic. Any treatment may only be performed if there is visual observation of the action of the unit.

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Safety features

The Endo-Arthroflator VET has the following safety features:

General safety functions

• **Monitoring of gas supply pressure:** A prolonged acoustic signal sounds (short triple peep with long pauses) if distension is activated at a gas supply pressure of < 10 bar or if the gas supply pressure drops below 10 bar during distension (red bottle symbol ③).

• **Negative alarm (in LAP mode):** If negative pressure is measured at the output to the patient, an acoustic signal sounds (quick double peep with long pauses, see picture on left). The warning continues until the pressure measured is positive again. This warning can also be taken as positive feedback. For example when the Veress needle is correctly positioned, the abdominal wall rises up and negative pressure arises in the abdomen. The negative alarm only sounds in LAP mode.

• **Software monitoring of pressure in the LAP mode:** If the patient pressure is ≥ 5 mmHg above the set setpoint an acoustic signal sounds after the set time delay (3 to 40 s, see service program 1, page 29) (long intermittent warning tone, see picture on left). If the gas flow is switched on, the ON/OFF solenoid valve is deactivated after 3 s. The gas flow is thus interrupted.

The acoustic alarm cannot be deactivated here. However, the acoustic signal can be set to sound after a relatively long time delay (e.g., 20 s), if this relatively long time delay has been set in the service program 1 (see page 29).

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



• **Software monitoring of pressure in the Arthro mode:** If the patient pressure is ≥ 5 mmHg above the set setpoint an acoustic signal (long intermittent warning tone - as in the LAP mode, see picture on left) sounds after a fixed time delay of 3 s. If the gas flow is switched on, the ON/OFF solenoid valve is also deactivated after 3 s. The gas flow is thus interrupted.

The acoustic alarm cannot be deactivated here. Remedy: See section above.

• **Hardware monitoring of pressure:** If the patient pressure is > 30 mmHg (LAP), >110 mmHg (ARTHRO) when the gas flow is switched on the ON/OFF solenoid valve is deactivated after a fixed time delay (approx. 5 s). The gas flow is interrupted and an acoustic signal sounds (continuous beep, see picture on left).


• **Mechanical monitoring of pressure:** If the patient pressure in the LAP mode is between 55 and 100 mmHg the blue overpressure safety valve opens. If the patient pressure in the ARTHRO mode is between 120 and 180 mmHg the uncolored overpressure safety valves opens. When the device is disconnected from the power supply, the LAP mode safety valve is connected.

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Error messages

- **Pressure sensor (E01):** One or both of the pressure sensors is/are defective when the device is switched on. An intermittent tone then sounds and the message E01 appears in the display.
- **Flow sensor (E02):** A self-test is performed for the flow sensor when the device is switched on. If this is not successful, an intermittent tone sounds and the message E02 appears in the display.
- **Checking the setpoints (E03):** If a false checksum is discovered in the non-volatile memory an intermittent tone sounds and the message E03 appears in the display. If one of the saved values is not within the permitted range, a default value (standard value) is taken.

Note: Make sure the device performs a self-check procedure at least once every day of operation.

All the light-up elements are firstly switched on. If no errors are found when checking the sensors and checksum this is acknowledged by 8 beeps.

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Unpacking the equipment

Carefully unpack the Endo-Arthroflator VET and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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
Basic equipment

- 1 Endo-Arthroflator VET 62432520
- 1 Insufflation line 8x12mm 20400143
- 1 Pack of sterile filters 031122-10 (mtp), disposable. Pack contains 10.
- 1 Power cord 400 A
- 1 Power cord 400 B
- 1 Instruction manual 96186008D
- 1 Cover for instruction manual
- 1 CO₂ supply tube, length as per order
- 1 Universal wrench 20400030

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Installing and connecting up the unit

Note: The Endo-Arthroflator VET including the accessories connected may be used only in medical facilities having electrical installations conforming to the applicable national, state, and local electrical codes.

Warning: The unit is not intended for use in hazardous zones. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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
The unit is equipped with a connector ⑭ for attaching a ground line.

The unit's ground line should be installed by a qualified electrician.

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Before connecting the unit to the electrical supply lines, verify that the line voltage stated on the unit's identification plate agrees with that of the electrical supply lines to be used.

Insert the power cord into the unit's power cord receptacle ⑯, pushing the cord's mating connector firmly home.

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Warning: The unit should only be connected to or disconnected from the power supply outside of areas subject to explosion hazards.

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Connection to a CO₂ gas bottle

Connect the gas output of the CO₂ bottle with high-pressure tube to the device's gas connection.

***Note:** The high-pressure tube is not supplied as standard, see Accessories p. 46 (available in 3 sizes).

Note: The gas connection is a standard American type. CO₂ bottles with PIN-Index or German-standard connections must be connected with an appropriate adapter tube.

CO₂ bottles with German/ISO connection

Connect the adapter tube to the outlet of the CO₂ bottle and screw tight using the universal wrench provided.

CO₂ bottles with PIN-Index connection

Place the connecting piece of the adapter tube on the CO₂ bottle and screw tight. To open the CO₂ bottle use the double open end tapped spanner provided.

Warning: Ensure that the CO₂ gas bottle is always upright when in use and that it is secured to prevent it from falling over (e.g. using a stand).

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Operating

When switched on, the device works exclusively in LAP mode. To switch the ARTHRO mode on, the device must be switched over (see first section on the ARTHRO mode on this page).

Switch on the device for the LAP mode

- Switch on the unit at its main power switch ①.
- After successful completion of the self-test 8 successive tones sound.
- The display flashes and must be confirmed. In order to confirm both setpoints (pressure and flow), it is sufficient to press the ± button ⑪ (pressure), the button ⑩ (flow) or the START/STOP button once.
- The device is now in LAP mode (intermittent gas flow type, i.e. with pumping pauses).
- The pressure setpoint can be set between 0 and 15 mmHg in increments of 1 mmHg.
- The flow setpoint can be set to the fixed programmed values 0 / 1.5 / 5 / 10 and 20 l/min.
- The flow setpoint can be increased in stages (e.g. from 0 to 1.5 and then from 1.5 to 5 etc) by pressing the button ⑩ (Flow). Once the maximum value (20) has been reached and the button ⑩ is pressed again, the device is set back to the minimum value (0).
- The gas volume display has a display range of 0 to 999 l (resolution 1 l). The display stops at 999. The display can be reset to 0 by pressing the Reset button ⑧.

Switch on the device for the Arthro mode

- Switch on the power switch ① and, at the same time, keep the + and - button for setting the pressure ⑪ held down for at least 3 seconds (display fields light up briefly). Once the term "ARTHRO" appears in the displays, the buttons + and - for setting the pressure can be released. The switch-over on the flashing value display is automatic.

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- The flashing display must now be confirmed. To this end, it is sufficient to press the ± buttons ⑪ (Pressure), the button ⑩ (Flow) or the START/STOP button once.
- The device is now in Arthro mode (continuous gas flow type, i.e. without pumping pauses).
- The Arthro mode is indicated by a patient pressure symbol ⑫ illuminated in red next to the + and - button for setting the pressure ⑪.
- The pressure setpoint can be set between 50 and 100 mmHg in increments of 5 mmHg
- The flow setpoint can be set to the fixed programmed values 0 / 1.5 / 3 / 5 and 10 l/min. The default value of 1.5 l/min is always set when the device is switched on.
- The flow setpoint can be increased in stages (e.g. from 1.5 to 3 and then from 3 to 5 etc) by pressing the button ⑩ (Flow). Once the maximum value (10) has been reached and the button ⑨ is pressed again, the device is set back to the minimum value (0).
- The gas volume display has a display range of 0 to 99 l (resolution 0.1 l). The display stops at 99. The display can be reset to 0 by pressing the Reset button ⑧. - In this mode an LED in the control panel is also illuminated

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Please pay attention to the status display for the gas input pressure ③.

Signal meanings:

- If the pressure values drop and a value of less than 10 bar is recorded, the display changes from green to red.
- If the pressure values increase and the 10 bar limit is exceeded, the display changes from red to green.

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Caution: To avoid contamination of the device as a result of the reverse flow of CO₂ or body fluid, a sterile CO₂ gas filter must be inserted between the insufflation tube and the device connection (Art. Nr. 031122-25 (mtp), for single use).

Place a sterile CO₂ gas filter on the insufflation connection ⑦. Fit insufflation hose to CO₂ gas filter.

Connect the other end of the insufflation tube to the desired instrument (e.g. VERESS needle, endoscope, trocar).

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Preparing for CO₂ insufflation

Enter the patient-line insufflation pressure desired using the “+” and “-” buttons ⑪.

The preselected patient pressure can be read from digital display ④. Blinking of the decimal point of digital display ④ indicates that the unit is in setpoint entry mode.

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Note: Brief, intermittent, depressions of the “+” and “-” buttons will gradually alter the pressure and flow rate setpoints. Holding these buttons depressed will rapidly alter these setpoints, allowing quickly reaching desired setpoint ranges.

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Use the button ⑩ to set the gas flow rate desired.

The preselected flow rate can be read from digital display ⑤. Blinking of the decimal point of digital display ⑤ indicates that the unit is in setpoint entry mode.

The fixed values are set in increments and ultimately the maximum value by repeatedly pressing the button ⑩. If this button is then pressed again, the display is reset to the minimum value 0.0 l/min.

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Note: Digital displays ④ and ⑤ will shift from indication of the values of the respective setpoint to indication of the current actual values of patient pressure flow rate when the insufflation START/STOP button ⑨ is pressed.

Note: If the insufflation START/STOP button ⑨ were to be pressed again now, insufflation would be started.

Pressing the button a second time would switch insufflation off.

Note: The setpoints can not be altered when insufflation is “On” (actual value display).

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
Note: The set pressure setpoints are saved when the device is switched off. The following applies for the flow: In LAP mode, the setpoint set prior to switching off is valid when the device is switched on again. In ARTHRO mode, the default value of 1.5 l/min is always set when the device is switched on.

When the device is switched on again, the values saved flash on the display and must be confirmed by pressing the ± button ⑪ for pressure, the button ⑩ for flow or the START/STOP button ⑨.

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Creation of the pneumoperitoneum

Start the insufflation process by pressing the Start/Stop (9) button. Open the connection valve to free the tube of air. The VERESS needle can now be introduced.

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Limitation of the flow by the small diameter of the VERESS needle is purely physical. This means, if the flow rate is set any higher, it will not be reached. In addition, the flow is controlled electronically in such a way that pressure conditions remain stable.

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The insufflation pressure, flow and gas volume can be read continuously on displays (4), (5) and (6).

If the patient pressure exceeds the preselected value by ≥ 5 mmHg, a long intermittent warning tone will sound. In the Arthro mode this occurs after a delay of precisely 3 seconds and in the LAP mode after a delay of 3 to 40 seconds (depending on the set value - see next page). A short intermittent warning tone will sound in the event of underpressure.

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Caution: If, during an arthroscopic intervention, both CO₂ and irrigation fluid are simultaneously connected, never open both stopcocks on the arthroscopy sheath at the same time, as this may cause a reverse flow of irrigation fluid into the unit and damage it.

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Service functions

In addition to the normal operating mode, the Endo-Arthroflator VET can also be started in a service mode in which certain additional functions can be selected. This mode is accessed by pressing the Reset button ⑧ when switching on the line voltage. Normal operation is no longer possible in the service mode (except in the service program 1). To exit this mode, the device must be switched off and on again.

Service program 1: Setting the time delay of the software alarm in LAP mode

After P01 has been displayed for a second, the programmed delay is then shown in seconds. The time delay for the software overpressure alarm can be set in increments of 1s from 3 to 40s with the flow setting button ⑩. If the maximum value (40) is set and the button ⑩ is pressed again, the minimum value (3) is displayed again. The Reset button ⑧ must be pressed to save the values. Once the values have been saved, the power is to be switched off and the device is to be switched on again at the power switch in normal mode.

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
Fuse replacement

Switch off the unit and remove the power plug from the electrical outlet.

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Release the line fuse holder ⑮ using a slim screwdriver (tip approx. 3 mm wide) or another suitable tool.

To this end, insert the tip of the screwdriver in the groove (see arrow) and, using the tip, push the clip below the groove down gently in order to release the line fuse holder from its fixation.

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Caution: Use only fuses that have the correct rating.

Insert new fuses of the appropriate rating.

	100...240 V~
Power fuse	2 x T1AL / 250 V

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Refit the power fuse holder ⑮.
Reconnect the power cord.
Test the unit for proper operation.

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Cleaning, disinfection and sterilization

Note: The instructions for “Cleaning, Sterilization and Care of KARL STORZ Instruments” must be followed. There the procedures for Cleaning, Disinfection and Sterilization are explained in detail.

Warning: Always pull out power plug before cleaning!
Caution: Avoid allowing fluids to enter the unit.

The exterior surfaces of the device can be cleaned with a disposable cloth moistened with disinfectant (no solutions containing alcohol).

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Inspection procedures for reusable silicone tubing sets before each sterilization

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Warning: A visual and functional check of the silicone tubing set must be carried out!

1. Visually inspect the silicone tubing set prior to each use to ensure that it is suitable for re-use. There should be no holes, cuts, dents or scratches in the silicone tubing. If the tubing set is damaged as described above, it should be replaced.
2. The integrity of the tubing set can be tested by attaching a syringe to one end of the tubing set and kinking the other end. Place the tubing in water. Force air into the tubing set using the syringe. If any bubbles appear in the water while forcing air through the tubing, the tubing set should be discarded.

Note: The leakage tester 13242 XL may also be used to check for leaks (see page 34).

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Cleaning, disinfection and sterilization of reusable tubing

Warning: Reusable tubing must be cleaned, disinfected and sterilized before use and subsequent reuse using validated preparation procedures.

Warning: Prior to the sterilization of reusable tubing, particularly in cases of manual preparation, it must be ensured that all chemical residue has been rinsed off or removed (see also preparation steps).

Warning: Disposable tubing sets, locks and hydrophobic bacteria filters must not be reused.

Note: For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the preparation of medical devices (see Appendix, pages 48, 49).

Preliminary cleaning of reusable tubing

Place the reusable tubing in a container filled with a cleaning solution (as per manufacturer's instructions) immediately after use to prevent contaminants from drying on the surface of the tubing set.

Caution: When preparing and using the solutions, follow the manufacturer's instructions paying close attention to proper dilution and working times. Prolonged immersion can lead to material changes.

Water Quality Requirements

Tap water can be used to prepare the cleaning and disinfectant solutions. When rinsing for the last time for example to remove any chemical residue, microbiologically pure/sterile water must be used.

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Preparation

Manual preparation of reusable tubing

Warning: For all cleaning and disinfection procedures on contaminated instruments, the local guidelines of the Employer's Liability Insurance Association or equivalent organizations must be observed.

1. Immerse reusable tubing in cleaning solution (as per manufacturer's instructions).
2. Using a syringe, rinse the tubing through underneath the (water) surface. Repeat the procedure several times if necessary. Remove contaminants with brushes, sponges, soft cloths or a cotton applicator.
Cleaning accessories are available from KARL STORZ (see catalog).
3. After cleaning, rinse the tubing and place it in a suitable container with disinfectant solution, making sure that all air bubbles have been removed. Ensure that the tubing is completely covered in the solution for the duration of exposure.
4. Cleaning brushes must be cleaned and disinfected immediately.
5. After disinfecting the tubing, remove it and rinse with microbiologically pure/sterile water. Ensure that all chemical residue in and on the tubing is removed.
6. Dry the tubing with sterile compressed air.
7. After preparation, inspect the tubes for cleanliness and damage. If soiling can still be detected, the preparation procedure must be repeated.

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Machine preparation of reusable tubing

To prepare using a machine, connect the reusable tubing to the retaining mechanisms in such a way that the instrument channels can be completely rinsed through.

Visible, heavy soiling must be removed in preliminary cleaning prior to machine preparation.

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Warning: Check the entire tubing set for leaks after preparation. As a result it is possible to detect a leak early, through which gas could escape.

To check for leaks use leak tester 13242 XL for example. For this purpose connect one end to the leak tester and block off the other end, with your finger for example. If the tubing set develops a leak, a pressure drop will be observed on the pressure gauge of the leakage tester.

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Sterilization

We recommend steam sterilizing the reusable tubing before use at 134 °C using the fractionated prevacuum procedure.

Warning: Successful sterilization is only possible on clean, disinfected surfaces.

Warning: The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and validated.

Warning: Any deviations from the recommended parameters for sterilization should be validated by the user.

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Steam Sterilization

KARL STORZ has validated the following sterilization method:

Fractionated prevacuum (steam sterilization)

Warning: Tubing sets may only be sterilized using the fractionated prevacuum procedure.

Trays should be positioned in the sterilizer so that there is adequate circulation and penetration of steam, air removal and condensate drainage.

Load the sterilizer according to its instruction manual.

At the completion of the steam sterilization cycle, all instruments should remain untouched until adequately cooled.

Warning: Burn hazard! Parts are hot after steam sterilization. Allow them to cool down.

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Fractionated prevacuum procedure

The fractionated prevacuum procedure consists of four phases.

Conditioning phase

The conditioning phase removes air from the chamber by pulling a vacuum up to four times. Steam is then blown in.

Sterilization phase

Sterilization takes place at a temperature of 134 °C (+3 °C) and an pressure (absolute) of 300 kpa (3 bar) over a period of 5 minutes (minimum time).

Evacuation phase

The steam is released.

Drying phase

Drying takes place under renewed application of a vacuum over a period of approx. 5–20 minutes.

Validated parameters for the fractionated prevacuum method:

Temperature:	134 °C (+3 °C)
Pressure (p _{absolute}):	300 kpa (3 bar)
Exposure time:	5 minutes for all instruments.

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
Maintenance

Regular maintenance can contribute to identifying potential problems before they become serious, enhancing the device's reliability and extending its useful operating life. Maintenance s can be obtained from your local representative or from the manufacturer.

Maintenance (functional and safety test) is recommended once a year, regardless of the accident or risk prevention regulations and/or testing intervals for medical instruments prescribed in different countries.

For detailed instructions please refer to the version of the Manual currently valid.

Note: To ensure safe operation of the unit, you should perform the following procedures every 12 months:

1. Check leakage current from chassis < 100 µA.
2. Check ground impedance < 0.1 ohms.

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Servicing and repair

Defective items of equipment must be d and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap. Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

For Endo-Arthroflator VET, individual repairs are necessary. Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparatur
Dr. Karl-Storz-Str. 34
78532 Tuttlingen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

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KARL STORZ GmbH & Co. KG
Abt. Reparatur
Dr. Karl-Storz-Str. 34
78532 Tuttlingen

Manufacturer's warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

Note: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparatur
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

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Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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Troubleshooting

Warning: Always unplug the unit before all maintenance work.

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Symptom:
Complete failure of the unit.

Possible causes:

- No power from the power line.
- Defective power fuse.
- Power cord is not properly connected to power cord receptacle.

Remedy:

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction manual. Use only fuses of the correct rating.
- Push power cord firmly into power cord receptacle.

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Symptom:
No gas flow or displayed input pressure is too low.

Possible causes:

- The connection to the device or gas supply is not correct.
- Insufflation not activated

Remedy:

- Check the connection to the device or gas supply. Have the gas supply checked.
- Press the insufflation START/STOP button.

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Symptom:
No pressure buildup.

Possible causes:

- Leakage in the system of gas lines.
- Failure in pressure regulation electronics.

Remedy:

- Check system of gas lines, particularly all fittings, for leakage and replace, if necessary.
- Send unit in for repair.

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Standard compliance (for 62432520)

According to EN60601-1

1990+A1:1993+A2:1995 Main section 3.
'Protection against risk of electric shock'

**According to EN 55011 / DIN EN 61326, EMC
requirements**

Limits and methods of measurement for noise
interference suppression of industrial, scientific
and high-frequency devices (ISM)

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This is not a medical device.

This device bears the CE mark.

A code number after the CE mark indicates the
responsible notified body.

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Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated as repairable by their respective manufacturers.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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

Plan of pneumatics for the Endo-Arthroflator VET 62 4325 20

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Type

- 1 = Cleaner
2 = Disinfectant
3 = Suitable for ultrasound baths
A = Not suitable for flexible endoscopes (according to manufacturers' instructions)
(A) = Only for the flexible endoscopes of the new generation (year of manufacture > 2005)

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Extract from the list of approved chemicals for the preparation of telescopes and endoscopic instrument sets

Manual preparation

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, prismatic light deflectors. For motors and their handles the appropriate manual provides detailed information.

The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to concentration and exposure times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Ecolab GmbH & Co. OHG	Sekusept Plus	1, 2, 3
	Sekusept Pulver Classic	1, 2, 3
	Ultraschall AF	1, 2, 3
Esteer Pharma GmbH	Ultraschall Aktiv	1, 2, 3
	Ultraschall Classic	1, 2, 3
Fresenius AG	Atid	1, 2
	Atid plus	2
Holifa Polska Sp. zo.o	Sporcid FF	1, 2
	Poisept Holifa	1, 2
Indeba Industria E Comercio LTDA	Letahdeido	2
Jose Collado S.A.	Darodor 4000 Liquido	1
	Darodor 9000	2
	Darodor Sinaldehyd 2000	2
Lysoform Dr. Rosemann GmbH	Aldasan 2000	1, 2
	Almyrol	1, 2
	Desoform	1, 2, 3
Medichem International	Lysoformin 3000	1, 2
	MedDis	2
	MediZyme	1
Merz Hygiene GmbH	Edisonite Super	1, 3
	Mucadont-LS	2, 3
	Mucadont-Zymaktiv	1, 3
METASYS Medizintechnik GmbH (formerly LONZA) Orochemie	Mucocit-T neu	1, 2, 3
Promagent AB	ID 50	2
	A 20	1, 3, A
Schumacher, Dr. GmbH	Wavacide	2
	Descoton Extra	1, 2
	Descoton forte	2
Schülke & Mayr GmbH	Perfektan Endo	2, 3
	Perfektan Neu	1, 2
	Perfektan TB	1, 2, 3, A
Steris	Gigasept	2, 3
	Gigasept AF forte	1, 2
	Gigasept FF	1, 2
Weigert, Dr. GmbH & Co.	Gigasept Instru AF	1, 2, 3, A
	Gigasept Med forte	1, 2
	Gigazyme	1, 3
Whiteley Industries PTY.Ltd	Lysotel V	2
	EnzyCare 2	1, 3
	neodisher LM 2	1, 3
Whiteley Industries PTY.Ltd	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher medizym	1, 3
Whiteley Industries PTY.Ltd	neodisher SeptoClean	1, A
	neodisher SeptoPreClean	1, 2
	Aidal Plus	2
Whiteley Industries PTY.Ltd	Matrix	1

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*** Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the life of the instrument set.

For flexible endoscopes, only immerse the examination sheath, not the housing/handle part.

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Type

- 1 = Cleaner
2 = Disinfectant
3 = Suitable for ultrasound baths
4 = Neutralizer
5 = Rinse aid
A = Not suitable for flexible endoscopes
(according to manufacturers' instructions)

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Machine preparation

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack.

For motors and their handles the appropriate manual provides detailed information.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Schülke & Mayr GmbH	Thermosept alka clean	1
	Thermosept DK	2
	Thermosept ED	2
	Thermosept ER	1
	Thermosept RKN-zym	1
Steris	Instru-Klenz	1, 3
	neodisher DuoClean	1, A
Weigert, Dr. GmbH & Co.	neodisher FA	1
	neodisher FA forte	1, 3
	neodisher Oxivario	1
	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher mediklar	5
	neodisher medizym	1
	neodisher SeptoClean	1, 2, A
	neodisher Septo DA	2
	neodisher Septo DN	2

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our department.

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: The Endo-Arthroflator VET should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the Endo-Arthroflator VET. The accessories and cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Table 200				
Accessories and cables which have been shown to comply with EN/IEC 60601-1-2:				
Type	Shielded	Length [m]	Ferrite	Used for
PA	No	> 3	No	Potential equalization
Power cord	No	3	No	Power supply

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Table 201

Guidance and manufacturer's declaration – electromagnetic emissions

Endo-Arthroflator VET model **62432520** is intended for use in the electromagnetic environment specified below.
The customer or user of the Endo-Arthroflator VET model **62432520** should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	Endo-Arthroflator VET model 62432520 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Endo-Arthroflator VET model 62432520 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 202			
Guidance and manufacturer's declaration – electromagnetic immunity			
Endo-Arthroflator VET model 62432520 is intended for use in the electromagnetic environment specified below. The customer or user of the Endo-Arthroflator VET model 62432520 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT * (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Complies <5% UT* (>95% dip in UT) for 0.5 cycle Complies 40% UT (60% dip in UT) for 5 cycles Complies 70% UT (30% dip in UT) for 25 cycles Complies <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

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Table 204 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
Endo-Arthroflator VET model 62432520 is intended for use in the electromagnetic environment specified below. The customer or user of the Endo-Arthroflator VET model 62432520 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to < 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms 3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Endo-Arthroflator VET model 62432520, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:</p> $d = [3.5/3]\sqrt{P} \quad 150 \text{ kHz to } < 80 \text{ MHz}$ $d = [3.5/3]\sqrt{P} \quad 80 \text{ MHz to } < 800 \text{ MHz}$ $d = [7/3]\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
<p>a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.</p> <p>b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

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Table 206

**Recommended separation distances between portable and mobile RF communications equipment and the
Endo-Arthroflator VET model 62432520**

The Endo-Arthroflator VET model 62432520 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

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To be completed by the supplier/importer:
Company stamp/signature:

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Tipo del aparato: Número de serie:

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Firma/Fecha:

To be filled out by instrument owner:
Return address/company stamp:

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KARL STORZ Endoscopy South Africa, (Pty) Ltd.
8th Floor Convention Tower
Cnr Heerengracht & Coen Steytler, Foreshore
Cape Town 8001, South Africa
P.O.Box 6061, Roggebaai Cape Town,
8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ ENDOSCOPY Kasachstan
Khodjanova 17
050060 Almaty, Kazakhstan
Phone/Fax: +7 72 72 49 43 63, +7 72 72 49 41 00
e-mail: KarlStorzKazakhstan@itte.kz

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51, Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.
#05 – 08 San Centre, 171 Chin Swee Road
Singapore 0316, Singapore
Phone: +65 6376-1066, Fax: +65 6376-1068
E-Mail: karlstorz@pacific.net.sg

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595
Phone: +61 (0)2 9490 6700, Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Beijing Representative Office
Room 610, China Life Tower
No. 6, Chaowai Street
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Chengdu Representative Office
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Guangzhou Representative Office
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Tel. No. +65 63761066, Fax. No. +65 63761068
Email : infoasia@karlstorz.com.sg
Email : serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
#06-07 e-Centre @ Redhill
Singapore 159471
Tel. No. +65 65325548, Fax No. +65 65323832
Email : infosp@karlstorz.com.sg
Email : servicesgp@karlstorz.com.sg

KARL STORZ Endoscopy Japan K. K.
Bois Hongo Building 6Fl, 3-42-5 Hongo
Bunkyo-ku, Tokyo 113-0033, Japan
Phone: +81 3 5802-3966, Fax: +81 3 5802-3988
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

KARL STORZ GmbH & Co. KG

Mittelstrasse 8
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de
Web: www.karlstorz.com